



## **els extra low shrinkage ® seal**

### **Technical Data**

|                                                  |                            |
|--------------------------------------------------|----------------------------|
| Consistency                                      | low viscous liquid         |
| Colour                                           | transparent / white opaque |
| Depth of curing (ISO 6874)                       | > 4.5 mm / > 1.5 mm        |
| Adhesion to bovine enamel (shear test)           | > 10 MPa                   |
| Flexural strength (DIN EN ISO 4049)              | > 100 MPa                  |
| Compressive strength<br>(analog DIN EN ISO 9917) | > 200 MPa                  |
| Conformity with DIN EN ISO 6874                  | fulfilled                  |
| Shelf life                                       | 4 years                    |